

myOSLchip Eraser

Tabletop light eraser for myOSLchip dosimeters

myOSLchip Eraser is a highly efficient tabletop light eraser compatible with myOSLchip dosimeters. myOSLchip Eraser erases up to 48 myOSLchip dosimeters in one load within 2 hours.

The myOSLchip dosimeters are opened and closed within the device to avoid exposure to light and to avoid mistakes of detectors and numbered covers.

RadPro

...Radiation Protection for the Radiation Professionals...



myOSLchip Eraser

Tabletop light eraser for myOSLchip dosimeters

Specifications

Capacity	Up to 48 myOSLchip dosimeter with BeO detectors
Sample carrier	Drawer carrier plate
Erasing LEDs	48 LEDs, blue
Erasing time	120 min (standard)
Voltage	90 – 250 V, 50/60 Hz
Fuse	2 A (T) micro fuse, delayed action fuse
Dimensions	250 mm x 200 mm x 410 mm
Transport case	530 mm x 260 mm x 330 mm
Weight	ca. 15.7 kg (net) ca. 22.7 kg (with transport case)
Certification	CE, ISO-9001

Related products

myOSLeraser



myOSLchip dosimeter



All information in this brochure is subject to technical changes without notice.

Rotunda Scientific Technologies® LLC - 3732 Fishcreek Road 913, Stow, OH 44224 USA - Info@RotundaSciTech.com - Phone: +1 (330) 906-3403 - www.rotundascitech.com